


6447992


Fig. 1

1.) ADDITION OF TYPE-SPECIFIC PEPTIDE MIXTURES TO THE WELLS



TEST DIRECTION (1->4)



TEST DIRECTION (4->1)



2.) INCUBATION 1h



3.) WASH 3 TIMES WITH 0.05 Tween 20/PBS


4.) ADDITION OF SERUM / INCUBATE FOR 1h



5.) TRANSFER SERUM TO NEXT / INCUBATE FOR 1h



6.) TRANSFER SERUM TO NEXT WELL / INCUBATE FOR 1h



7.) TRANSFER SERUM TO NEXT WELL / INCUBATE FOR 1h



8.) WASH 3 TIMES WITH 0.05 % Tween 20 / PBS

9.) ADD SECOND ANTIBODY

10.) ADD ABTS SOLUTION: COLOUR REACTION

11.) MEASURE OD